

## Product Datasheet

### YW1159 (orb611994)

<b>Description</b>	YW1159 is an inhibitor of Wnt signaling with an IC50 value of 1.2 nM in a reporter assay.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{19}H_{14}FN_5O$
<b>MW</b>	347.3
<b>CAS Number</b>	[2131223-69-1]
<b>SMILES</b>	<chem>O=C(C1=C(C)N(C2=C(F)C=CC=C2)N=N1)NC3=CC=C(C=CC=C4)C4=N3</chem>
<b>Expiration Date</b>	12 months from date of receipt.