

Product Datasheet

Mouse PSMA / FOLH1 Protein (orb612091)

Description	Mouse PSMA, His Tag is expressed from human 293 cells (HEK293). It contains AA Ile 44 - Ala 752 (Acc...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PSMA
Preservatives	MES and NaCl
Form/Appearance	Powder
Storage	-20°C
Tag	N-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 81.5 kDa. The protein migrates as 90-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_058050.3
Purity	95%
Source	Mouse PSMA, His Tag (orb612091) is expressed from human 293 cells (HEK293). It contains AA Ile 44 - Ala 752 (Accession # O35409-1).
MW	81.5 kDa
Uniprot ID	O35409
NCBI	NP_058050.3
Expiration Date	6 months from date of receipt.