

Product Datasheet

Mouse TSLP Protein (orb612095)

Description	Mouse TSLP, His Tag is expressed from human 293 cells (HEK293). It contains AA Tyr 20 - Glu 140 (Acc...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TSLP
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 15.9 kDa. The protein migrates as 23-27 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_067342.1
Purity	90%
Source	Mouse TSLP, His Tag (orb612095) is expressed from human 293 cells (HEK293). It contains AA Tyr 20 - Glu 140 (Accession # Q9JIE6-1).
MW	15.9 kDa
Uniprot ID	Q9JIE6
NCBI	NP_067342.1
Expiration Date	6 months from date of receipt.