

Product Datasheet

Human CRTAM Protein (orb612121)

Description	Human CRTAM, His Tag is expressed from human 293 cells (HEK293). It contains AA Ser 18 - Gly 287 (Ac...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CRTAM
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 31.9 kDa. The protein migrates as 45-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_062550.2
Purity	90%
Source	Human CRTAM, His Tag (orb612121) is expressed from human 293 cells (HEK293). It contains AA Ser 18 - Gly 287 (Accession # O95727-1).
MW	31.9 kDa
Uniprot ID	O95727
NCBI	NP_062550.2
Expiration Date	6 months from date of receipt.