

Product Datasheet

Mouse Glypican 1 / GPC1 Protein (orb612132)

Description	Mouse Glypican 1, His Tag is expressed from human 293 cells (HEK293). It contains AA Asp 24 - Ser 52...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Glypican 1
Preservatives	50 mM Tris, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 57.5 kDa. The protein migrates as 60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_057905.1
Purity	90%
Source	Mouse Glypican 1, His Tag (orb612132) is expressed from human 293 cells (HEK293). It contains AA Asp 24 - Ser 529 (Accession # Q9QZF2-1).
MW	57.5 kDa
Uniprot ID	Q9QZF2
NCBI	NP_057905.1
Expiration Date	6 months from date of receipt.