

---

## Product Datasheet

### Human Glypican 1 / GPC1 Protein (orb612137)

**Description**

Human Glypican 1, Fc Tag is expressed

contains AA Asp 24 - Ser 530...

**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

Glypican 1

**Preservatives**Tris with Glycine, Arginine and NaCl,  
pH7.5**Form/Appearance**

Powder

**Storage**

-20°C

**Tag**

C-Fc

**Note**

For research use only

**Application notes**

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 82.5 kDa. The protein migrates as 90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Protein Sequence**

NP\_002072.2

**Purity**

90%

**Source**

Human Glypican 1, Fc Tag (orb612137) is expressed from human 293 cells (HEK293). It contains AA Asp 24 - Ser 530 (Accession # P35052-1).

**MW**

82.5 kDa

**Uniprot ID**[P35052](#)**NCBI**[NP\\_002072.2](#)**Expiration Date**

6 months from date of receipt.