

Product Datasheet

Human PVRIG Protein (orb612146)

Description	Human PVRIG, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Thr 41 - Asp 171 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PVRIG
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 40.3 kDa. The protein migrates as 45-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_076975.2
Purity	95%
Source	Human PVRIG, Fc Tag (orb612146) is expressed from human 293 cells (HEK293). It contains AA Thr 41 - Asp 171 (Accession # Q6DKI7-1).
MW	40.3 kDa
Uniprot ID	Q6DKI7
NCBI	NP_076975.2
Expiration Date	6 months from date of receipt.