

Product Datasheet

Human GITR Ligand / TNFSF18 Protein (orb612155)

Description

Human GITR Ligand, His Tag is

(HEK293). It contains AA Gln 50 - Ser 177.

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

GITR Ligand

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

N-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 20.5 kDa. The protein migrates as 23 kDa and 26 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

Gln 50 - Ser 177

Purity

90%

Source

Human GITR Ligand, His Tag (orb612155) is expressed from human 293 cells (HEK293). It contains AA Gln 50 - Ser 177 (Accession # AAH69319.1).

MW

20.5 kDa

Uniprot ID

AAH69319.1

Expiration Date

6 months from date of receipt.