



Product Datasheet

FLAG Tag (CT) antibody (orb612171)

Description	FLAG Tag (CT) antibody
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Flag Tag peptide
Target	FLAG Tag (CT)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB=1:5000-20000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: DDDD-K-Tagged Fusion Protein



Western blot: Lane 1: DDDD-K-Tagged Fusion Protein