
Product Datasheet

Alpha Synuclein antibody (pSer129) (orb612725)

Description

Rabbit monoclonal antibody against

conjugated to Biotin

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Biotin

Tested Applications

ELISA, ICC, IHC, WB

ImmunogenHuman alpha synuclein AA 124-134:
AYEMP-pS-EEGYQ-Cys**Target**

Alpha Synuclein pSer129

PreservativesPBS pH 7.4, 50% glycerol, 0.09%
Sodium azide *Storage buffer may
change when conjugated**Concentration**

1 mg/ml

StorageConjugated antibodies should be stored
at 4°C**Note**

For research use only

Application notesA 1:500 dilution of SMC-600 was
sufficient for detection of Alpha
Synuclein pSer129 in 10 µg of Mouse
Brain by ECL immunoblot analysis using
Goat Anti-Rabbit IgG:HRP as the
secondary antibody.**Isotype**

IgG

Clonality

Monoclonal

Uniprot ID[P37840](#)**NCBI**[NP_000336.1](#)**Entrez**[6622](#)**Dilution Range**

WB (1:500)

Expiration Date

12 months from date of receipt.