

Product Datasheet

mCherry Antibody (orb613681)

Description	Our Anti-mCherry primary antibody from Biorbyt is mouse monoclonal. It detects mCherry and is Protei...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Recombinant full length protein
Target	mCherry
Preservatives	PBS + 50% glycerol and 5 mM NaN ₃
Form/Appearance	Protein G Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	28
Dilution Range	WB: 1:1000, IHC: 1:500, ICC: 1:500
Expiration Date	12 months from date of receipt.