

Product Datasheet

CD3E antibody (orb613774)

Description	Rabbit monoclonal antibody to CD3E
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human peripheral E+ cells.
Target	CD3 epsilon
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P07766
Expiration Date	12 months from date of receipt.