

Product Datasheet

Anti-CD169 [3D6.112 (recombinant version)] (orb613783)

Description	Rabbit monoclonal antibody to CD169_Sialoadhesin
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Spleen cells from AO rats immunised with thioglycollate-elicited murine peritoneal macrophages (TPM) that had been cultured in mouse serum to induce sheep erythrocyte receptor (SER).
Target	CD169
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rat IgG
Clonality	Monoclonal
Clone Number	3D6.112 (recombinant version)
Purity	Purified
Uniprot ID	Q62230

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Expiration Date

12 months from date of receipt.



Anti-CD169 staining of murine macrophages Murine bone marrow-derived macrophages (BMDMs) were stained with anti-CD169 antibody orb613784 (as well as isotype control) and bound antibody detected using goat IgG anti-rat IgG (H&L-chain) polyclonal antibody directly conjugated to Alexa Fluor® 647 (AF647) commercially available from a competitor. Compared to the isotype control and secondary-only control, a CD169-positive population is apparent.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)