

Product Datasheet

CD169_Sialoadhesin antibody (orb613784)

Description	Rat monoclonal antibody to CD169_Sialoadhesin
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Spleen cells from AO rats immunised with thioglycollate-elicited murine peritoneal macrophages (TPM) that had been cultured in mouse serum to induce sheep erythrocyte receptor (SER).
Target	CD169
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q62230
Expiration Date	12 months from date of receipt.