

## Product Datasheet

### TRP-1 antibody (orb613791)

<b>Description</b>	Mouse monoclonal antibody to TRP-1
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC, IF, IHC-Fr, IP, WB
<b>Immunogen</b>	70-75 kDa pigmentation-associated glycoprotein in human melanoma cell lines.
<b>Target</b>	TRP-1, gp75
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	This antibody binds to human and murine TRP-1 and had been under evaluation for the treatment of tumours.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Uniprot ID</b>	<b>P17643</b>
<b>Expiration Date</b>	12 months from date of receipt.