

## Product Datasheet

### IFN gamma antibody (orb613796)

<b>Description</b>	Rabbit monoclonal antibody to IFN gamma
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Monkey
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	Blocking, ELISA, FC, WB
<b>Immunogen</b>	FPLC-purified natural human IFN- $\gamma$ . Refer to the original article [PMID: 3114127]
<b>Target</b>	IFN gamma
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	This antibody can be used to detect both native and recombinant intereferon gamma (IFN $\gamma$ ). It is reported only to bind to the functional, dimeric form of IFN $\gamma$ .
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Uniprot ID</b>	<b>P01579</b>
<b>Expiration Date</b>	12 months from date of receipt.