

Product Datasheet

CDCrel-1 antibody (orb613816)

Description	Mouse monoclonal antibody to CDCrel-1
Species/Host	Mouse
Reactivity	Human, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Immunoprecipitate from human brain.
Target	CDCrel-1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q99719
Expiration Date	12 months from date of receipt.