
Product Datasheet

IgE antibody (orb613864)

Evolution Bioproducts.
Description

Mouse monoclonal antibody to IgE

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA
Immunogen	Polyclonal human IgE.
Target	IgE
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	<p>The antibody binds specifically to IgE, an immunoglobulin isotype secreted by upon stimulation of B cells by interleukins IL-4 and IL-13. IgE provides immunity to parasites such as helminths and Plasmodium falciparum. The immunoglobulin also plays a central role in the development of Type I hypersensitivity reactions and allergic diseases such as asthma, sinusitis and food allergies. These allergic responses are triggered by the binding of IgE to the high affinity receptor FcεRI, expressed on the surface of basophils and mast cells. Binding of IgE causes cross-linking of the receptor, inducing degranulation and release of inflammatory mediators, particularly of histamine. The monoclonal antibody TES-C21 binds to secreted IgE and to IgE expressed on the surface of IgE-producing B cells, thereby inhibiting IgE binding to FcεRI and the release of histamine from mast cells. The monoclonal antibody can have diagnostic applications, such as for identifying and quantifying IgE-expressing B cells and IgE levels in serum samples, and therapeutic applications in the treatment of allergies. Furthermore, the antibody can be conjugated to cytotoxic drugs and deliver them specifically to IgE-producing B cells.</p>
Isotype	IgG1
Clonality	Monoclonal