

Product Datasheet

Bet v 1 antibody (orb613896)

Description	Mouse monoclonal antibody to Bet v 1
Species/Host	Mouse
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Bet v 1-A, a major allergen of the European white birch.
Target	Bet v 1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody binds to Bet v 1 via an epitope commonly used by allergogenic IgE. Fab versions of this antibody have been used in co-crystallization studies of Bet v 1.
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P15494
Expiration Date	12 months from date of receipt.