

Product Datasheet

Anti-PD-L2 [TY25] (orb613949)


Description	Rat monoclonal antibody to PD-L2
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC, IP, WB
Immunogen	This antibody was raised by immunising SD rats with B7-DC-transfected RAW264.7 cells and selected on B7-DC-transfected NRK.
Target	PD-L2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Clone Number	TY25
Purity	Purified
Uniprot ID	Q9WUL5
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Binding curve of anti-PDL2 antibody TY25 (orb613948) to mouse PDL2-Fc fusion protein ELISA Plate coated with mouse PDL2-Fc fusion protein at a concentration of 5 µg/ml. A 3-fold serial dilution from 10000 to 0.1 ng/ml was performed using orb613948. For detection, a 1:4000 dilution of HRP-labelled goat anti-rabbit antibody was used.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)