

Product Datasheet

PDHA1/2 (Phospho-Ser293/291) Antibody (orb614539)

Description	PDHA1/2 (Phospho-Ser293/291) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P08559
Expiration Date	12 months from date of receipt.