

Product Datasheet

c-SRC/FYN/c-YES (phospho-Tyr419/420/426) Antibody (orb615020)

Description	c-SRC/FYN/c-YES (phospho-Tyr419/420/426) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P12931
Expiration Date	12 months from date of receipt.