

Product Datasheet

HSP20 (phospho-Ser16) Antibody (orb615035)

Description	HSP20 (phospho-Ser16) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O14558
Expiration Date	12 months from date of receipt.