



Product Datasheet

P-V-Myb+C-Myb(Phospho-S11) Antibody (orb615064)


Description	P-V-Myb+C-Myb(Phospho-S11) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P10242
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of human ut...



Immunohistochemical analysis of human br...



Immunofluorescence analysis of A9 cells...