

Product Datasheet

JAK2(Phospho-Y1007+Y1008) Antibody (orb615070)

Description	JAK2(Phospho-Y1007+Y1008) Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	O60674
Expiration Date	12 months from date of receipt.