



Product Datasheet

p95/NBS1 (Phospho-S343) Antibody (orb615072)

Description	p95/NBS1 (Phospho-S343) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O60934
Expiration Date	12 months from date of receipt.



ICC staining Phospho-p95/NBS1(S343) in P...



Immunohistochemical analysis of paraffin...