



## Product Datasheet

### GSK3 beta(Phospho-Ser 9) Antibody (orb615074)


<b>Description</b>	GSK3 beta(Phospho-Ser 9) Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IHC, WB
<b>Storage</b>	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P49841</b>
<b>Expiration Date</b>	12 months from date of receipt.



ICC analysis of hela cells (green). The ...



ICC analysis of MCF-7 cells (green). The...



Immunohistochemical analysis of paraffin...