



Product Datasheet

Cdc6 (Phospho-S54) Antibody (orb615132)


Description	Cdc6 (Phospho-S54) Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q99741
Dilution Range	WB: 1:1,000-1:2,000, IF/ICC: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Phospho-Cdc6(S54) in HeLa cells.



Immunofluorescence analysis of HeLa cells.



Immunofluorescence analysis of HepG2 cells.