

Product Datasheet

Rac1+Cdc42(Phospho-Ser71) Antibody (orb615143)

Description	Rac1+Cdc42(Phospho-Ser71) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P60953
Expiration Date	12 months from date of receipt.