



Product Datasheet

NEXN Polyclonal Antibody (orb615847)

Description	NEXN Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human NEXN
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
Purity	Affinity purification
MW	81kDa
Uniprot ID	Q0ZGT2
Dilution Range	WB 1:500 - 1:2000; IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of U-2 OS ce...