

Product Datasheet

Acetyl and phospho-Histone H3 (Ac-K9/p-Ser10) antibody (orb6162)

Description

Acetyl and phospho-Histone H3 (Ac-K9/p-Ser10)

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Histone H3 around the site of Acetyl-Lys9 and phospho-Ser10 (AR(Ac-K)(p-S)TG)
Target	Histone H3 (acetyl K9, phospho S10)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	15 kDa
CAS Number	Histone H3.1
Uniprot ID	P68431
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Pancreas (Rat); Lane ...