

Product Datasheet

Acetyl-Histone H2B-K12 pAb Antibody (orb616709)

Description	Acetyl-Histone H2B-K12 pAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q16778
Expiration Date	12 months from date of receipt.