


Product Datasheet

ICAM1 (Phospho-Tyr512) antibody (orb6189)

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to ICAM1 (Phospho-Tyr512) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| Dilution Range | WB:1:500-1:2000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of HeLa whole cell...