

Product Datasheet

Mouse anti Rabbit IgG heavy chain(Fc)-HRP Antibody (orb620951)

Description	Mouse anti Rabbit IgG heavy chain(Fc)-HRP Antibody
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Expiration Date	12 months from date of receipt.