



Product Datasheet

IRX2 Antibody (orb621037)

Description	IRX2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9BZ11
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cell lysat...



Immunofluorescent staining of HeLa cells...