


## Product Datasheet

### SNAP-25 Antibody (orb621147)

<b>Description</b>	SNAP-25 Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IP, WB
<b>Storage</b>	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P60880</b>
<b>Dilution Range</b>	WB: 1:1000, IF/ICC: 1:100, IP: 1:10-100
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis on SH-SY-5Y cells ...



Immunofluorescent staining of SNAP25 in ...