




Product Datasheet

Myoglobin Antibody (orb621770)

Description	Myoglobin Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P02144
Dilution Range	WB: 1:1,000-5,000, IHC-P: 1:50-1:200, IF/ICC: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Myoglobin on hu...



Immunohistochemical analysis of paraffin...