


Product Datasheet

IRAK1 (phospho-Ser376) antibody (orb6221)

Description	IRAK1 (phospho-Ser376) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human IRAK1 around the phosphorylation site of Ser376 (QS(p-S)TV)
Target	Phospho-IRAK1 (Ser376)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	78 kDa
Uniprot ID	Q62406
Dilution Range	IHC-P=1:100-500, ICC=1:100, Flow-Cyt=1ug/test, WB=1:500-2000, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 lysates us...