



Product Datasheet

CD19 Antibody (orb622142)

Description	CD19 is a transmembrane glycoprotein that contains two extracellular immunoglobulin-like domains. CD...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	The immunogen is Synthetic peptide within Human CD19 aa 384-430. The mAb will bind to the intracellular region
Form/Appearance	1*TBS (pH7.4), 1%rAlbumin, 40%Glycerol. Preservative: 0.05% Sodium Azide
Storage	Store at -20°C
Note	For research use only
Clonality	Recombinant
Purity	ProA affinity purified
MW	75-100 kDa
Uniprot ID	P15391
Dilution Range	WB: 1:1,000-5,000, FC: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of CD19 on different samples.



000 dilution. Positive control: Lane 1: ...