



Product Datasheet

MMP9 Antibody (orb622328)


Description	MMP9 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P14780
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of A431 cells u...



ICC staining HeLa cells (green). The nuc...



ICC staining SHG-44 cells (green). The n...