

## Product Datasheet

### CD137 antibody (APC) (orb623063)

<b>Description</b>	Mouse monoclonal antibody against CD137 conjugated to APC
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	APC
<b>Tested Applications</b>	FC
<b>Immunogen</b>	Recombinant human CD137 ectodomain
<b>Target</b>	CD137
<b>Preservatives</b>	Purified IgG conjugated to Allophycocyanin (APC) - liquid
<b>Form/Appearance</b>	Purified IgG conjugated to Allophycocyanin (APC) - liquid
<b>Storage</b>	+4°C DO NOT FREEZE
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">Q07011</a>
<b>Expiration Date</b>	12 months from date of receipt.