

Product Datasheet

Goat anti Rat IgG antibody (orb623142)

Description	Goat polyclonal anti Rat IgG
Species/Host	Goat
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, WB
Immunogen	Normal rat IgG
Target	IgG
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.