

Product Datasheet

Goat anti Rat IgG antibody (Biotin) (orb623143)

| | |
|----------------------------|---|
| Description | Goat polyclonal anti Rat IgG conjugated to Biotin |
| Species/Host | Goat |
| Reactivity | Rat |
| Conjugation | Biotin |
| Tested Applications | ELISA, IF, WB |
| Immunogen | Normal rat IgG |
| Target | IgG |
| Preservatives | Purified IgG conjugated to Biotin - liquid |
| Form/Appearance | Purified IgG conjugated to Biotin - liquid |
| Concentration | 0.5 mg/ml |
| Storage | +4°C, -20°C if preferred |
| Note | For research use only |
| Isotype | Polyclonal IgG |
| Clonality | Polyclonal |
| Expiration Date | 12 months from date of receipt. |