

Product Datasheet

Goat anti Rat IgG antibody (FITC) (orb623144)

Description	Goat polyclonal anti Rat IgG conjugated to FITC
Species/Host	Goat
Reactivity	Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	Normal rat IgG
Target	IgG
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.