

Product Datasheet

Goat anti Rat IgG antibody (HRP) (orb623145)

Description	Goat polyclonal anti Rat IgG conjugated to HRP
Species/Host	Goat
Reactivity	Rat
Conjugation	HRP
Tested Applications	ELISA, IH-Fr, WB
Immunogen	Normal rat IgG
Target	IgG
Preservatives	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.