

Product Datasheet

Goat anti Rabbit IgG/IgM antibody (Biotin) (orb623147)

Description	Goat polyclonal anti Rabbit IgG/IgM conjugated to Biotin
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Biotin
Tested Applications	ELISA, IH-Fr, IHC-P, WB
Immunogen	Normal rabbit IgG + IgM.
Target	IgG IgM
Preservatives	Purified IgG conjugated to Biotin - liquid
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	IgG concentration 0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Uniprot ID	P01780, P04221, P03988
Expiration Date	12 months from date of receipt.