

Product Datasheet

Goat anti Rabbit IgG/IgM antibody (HRP) (orb623148)

Description	Goat polyclonal anti Rabbit IgG/IgM conjugated to HRP
Species/Host	Goat
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, IH-Fr, IHC-P
Immunogen	Normal rabbit IgG + IgM.
Target	IgG IgM
Preservatives	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Uniprot ID	P01780, P04221, P03988
Expiration Date	12 months from date of receipt.