

## Product Datasheet

### Goat anti Rabbit IgG/IgM antibody (Texas Red) (orb623150)

<b>Description</b>	Goat polyclonal anti Rabbit IgG/IgM conjugated to Texas Red
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	Texas Red
<b>Tested Applications</b>	FC, IHC-Fr, IHC-P
<b>Immunogen</b>	Normal rabbit IgG + IgM.
<b>Target</b>	IgG IgM
<b>Preservatives</b>	Purified IgG conjugated to Texas Red - liquid
<b>Form/Appearance</b>	Purified IgG conjugated to Texas Red - liquid
<b>Storage</b>	+4°C, -20°C if preferred
<b>Note</b>	For research use only
<b>Isotype</b>	Polyclonal IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P01780</a> , <a href="#">P04221</a> , <a href="#">P03988</a>
<b>Expiration Date</b>	12 months from date of receipt.