

Product Datasheet

Recombinant Rat IL-17A (orb623210)

Description	IL-17A is a member of the IL-17 family of cytokines, whose members are involved in numerous immune regulatory functions. IL-17 induces the production of many other cytokines, chemokines, and prostaglandins. Rat IL-17A Recombinant Protein is purified interleukin-17A produced in yeast.
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Application notes	The Mouse M-CSF protein can be used in cell culture, as a M-CSF ELISA Standard, and as a Western Blot Control.
Protein Sequence	AVLIQSSVC PNAEANNFLQ NVKVNLIKVLN SLSSKASSRR PSDYLNIRSTS PWTLSRNEDP DRYPSVIWEA QCRHQRCVNA EGKLDHHMNS VLIQQEILVL KREPEKCPFT FRVEKMLVGV GCTCVSSIVR HAS (133)
Source	Yeast
MW	15.0 kDa
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)