

Product Datasheet

Cytomegalo Virus Pp28 (UL99) Antibody (orb623311)

Description	Cytomegalo Virus Pp28 (UL99) Antibody
Reactivity	Human
Tested Applications	ELISA, WB
Immunogen	Recombinant protein contains the CMV Pp28 (UL99) immunodominant regions, 130-160 amino acids
Form/Appearance	Liquid
Storage	Store at room temperature for up to 3 weeks. Store desiccated at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.